

Image Acquisition - Sensors and Sources

William W. Moses

Lawrence Berkeley National Laboratory

March 17, 2003


Outline:

- **Introduction**
- **Benefits of Improved Sensors**
- **Recommendations**

Existing Sensors Limit Performance!!!

Example: PET

ing of Photon
detectors



- Radiotracer decays by emitting a positron (β^+).
- 511 keV photon pairs detected via time coincidence.
- Positron localized to line (defined by detector pair).
- Image reconstructed using computed tomography.

- Perfect Sensor \Rightarrow Time-of-Flight Localizes Positron to Point
 - PET becomes PI (Positron Imaging)?


Energy per Quantum

| <u>Imaging Modality</u> | <u>Energy / Quantum</u> |
|-------------------------|--------------------------------|
| PET | $5 \times 10^5 \text{ eV}$ |
| SPECT | $1 \times 10^5 \text{ eV}$ |
| X-Ray | $5 \times 10^4 \text{ eV}$ |
| Optical | $2 \times 10^0 \text{ eV}$ |
| Ultrasound | $2 \times 10^{-2} \text{ eV}$ |
| MRI | $1 \times 10^{-13} \text{ eV}$ |


Wide Range of Energies \Rightarrow Wide Range of Sensors

Sensors for Medical Imaging

Scintillation Detector




Solid State (or Direct) Detector



**Sensors Are Optimized for Each Modality,
But General Concepts are Identical**

Typical Sensor for PET

Array of Scintillator Crystals
($\text{Bi}_4\text{Ge}_3\text{O}_{12}$ or $\text{Lu}_2\text{SiO}_5:\text{Ce}$)



Photomultiplier Tubes

PET Needs

Improve Timing Resolution

- Reduce Random Event Background (Up to 50% of Total Events)
- Use Time-of-Flight to Reduce Noise Variance (by 5x-10x)

Improve Energy Resolution

- Scattered Events Often Outnumber True Events

Improve Spatial Resolution

- Isotropic Resolution Throughout the Field of View
- Small Animal Imaging

SNR is Main Limitation in Clinical PET

PET: Potential Sensor Improvements




New Scintillators

- Improved Timing Resolution
- Improved Energy Resolution

Avalanche Photodiodes

- Smaller Pixels improves Spatial Resolution
- Higher Quantum Efficiency (80% vs. 25%), can improve Energy Resolution
- Better Stability



Solid State Detectors Not Useful for PET

- Too Slow (speed of sound vs. speed of light)
- High Cost per Unit Volume

New Sensors Can Improve Spatial Resolution and Significantly Reduce Background

Typical Sensor for SPECT



**Photomultiplier
Tubes**

**Scintillator
Crystals
(NaI:TI)**

SPECT: Needs

Improve Energy Resolution

- Reduce Scatter Background
- Over 35% of SPECT Events are Scatter

Improve Stability


- Maintain Calibration

Improve Resolution / Efficiency Tradeoff

- Set by Collimator
- Major Improvements Difficult


Energy Resolution is Main Limitation in SPECT

SPECT: Potential Sensor Improvements





New Scintillators

- Improved Energy Resolution




Solid State Detectors

- *Much* Improved Energy Resolution
- Improved Stability



Best Theoretical Energy Resolution with Solid State

Typical Sensor for X-Ray CT



X-Ray CT: Needs

Reduced Dose

Improved Spatial Resolution

- 0.5 mm to 0.2 mm
- Difficult Without Increasing Dose or Increasing Noise

Reduced Imaging Time

- Factor of 2

Scintillator & Photodetector: Little Improvement Possible



- Scintillator: High Luminous Efficiency & QE, Low Afterglow
- Photodetector: High QE, Very Stable

Main Limitation in CT is Dose

X-Ray CT: Potential Sensor Improvements

Sensor: Larger 2-D Arrays

- High Channel Density
(Thermal Management,
Connection Density)
- Maintain Mechanical Tolerances
- Want Thinner Reflector




Source: Higher Power

- Same Dose in Shorter Time
- Monochromatic Source?



Significant Reduction in Acquisition Time

Typical Sensor for Planar X-Ray


Film / Screen



Digital




Planar X-Ray: Needs



Scintillator & Photodetector: Little Improvement Possible

- High Luminous Efficiency
- High Detective Quantum Efficiency ($>70\%$)




Some Improvement Possible:

- More Portable, Durable, Lighter
- Improved Falloff of DQE with Spatial Frequency
- Better Scatter Rejection / Dose Tradeoff?
 - Slot Scan?
 - More Efficient Scatter Grids?


Digital X-Ray Sensors Close to Ideal

Typical Sensor for Ultrasound



Piezoelectric Ceramic

Ultrasound Transducer Needs



Real-Time 3-D Imaging

- 2-D Arrays

Greater Bandwidth

- Same Transducer Used for Many Frequencies

Reduce Reflection at Tissue Interface

- Match Transducer & Tissue Impedance

Higher Frequency

- Higher Resolution \Rightarrow Small Animal Imaging

Room for Improvement in Many Aspects

Recommendation 1: New Materials

Material Discovery




- Scintillators
- Solid State Detectors
- Ultrasound Transducers
- *Very Little Current Effort*

Material Development

- Difficult & Expensive
(Minimum 10 years @ \$1M / year for each Material)
- Done by Commercial Sector (if Patent Protected)
- Done by Government (if Not Patent Protected)


rtually All Modalities Need New Sensor Material

Recommendation 2: Electronics & Packaging



All Modalities Want:

- Larger Area Covered
- Smaller Pixels
- ⇒ Many Pixels (25 million)




Implications:


- ⇒ High Power Consumption & Density
- ⇒ Large Number & Density of Connections
- ⇒ Modular Design
- ⇒ Dead Area Between Modules
- ⇒ Scintillator Arrays w/ Small Pixels



Packaging Improvements Would Help All Modalities



Thanks To:



Tom Budinger, *LBNL*
Mike Casey, *CPS Innovations*
Simon Cherry, *UC Davis*
Neal Clinthorne, *Univ. of Michigan*
Jim Colsher, *GE Medical Systems*
Steve Derenzo, *LBNL*
Ruvyn Deych, *Analogic*
Kathy Ferrara, *UC Davis*
Stuart Foster, *Univ. of Toronto*
Daniel Gagnon, *Picker Imaging*

Paul Granfors, *GE Medical Systems*
Grant Gullberg, *LBNL*
Anne Hall, *GE Medical Systems*
Floris Jansen, *GE Medical Systems*
Ron Keyser, *EG&G Ortec*
Jim LeBlanc, *GE Medical Systems*
Chuck Melcher, *CPS Innovations*
Kanai Shah, *RMD, Inc.*
Ling Shao, *Philips Medical Systems*

