

Stabilized Adaptive Optics Imaging for Laser Microsurgery

NIH NIBIB
Image Guided Interventions Workshop

D. X. Hammer
Physical Science Inc.
Andover MA

Support by EB003111

13-14 May 2004




Program Objectives

- **Overall Program Goal**
 - Design, construct, and test a prototype clinical imaging and irradiation system, based upon **adaptive optics, retinal tracking, and ultrashort laser pulse-induced retinal damage mechanisms.**
- **Specific Program Objectives**
 - Design and construct an adaptive optics scanning laser ophthalmoscope (AOSLO).
 - Integrate a retinal tracker into the AOSLO (TAOSLO).
 - Characterize the TAOSLO imaging and tracking performance in human subject tests at PSI.
 - Deliver and setup the TAOSLO for ultrashort pulse delivery at the Optical Radiation Branch (AFRL/HEDO, Brooks AFB TX) and test in an animal model.



Motivation and Clinical Significance


- Advances in diagnostic laser technologies have far outpaced advances in laser therapies (except for PDT)
- Retinal diseases
 - Age-related macular degeneration (AMD)
 - Diabetic retinopathy (DR)
 - Retinitis pigmentosa (RP)
 - Intraocular melanoma (IM)
- Potential therapies with the proposed system - Early detection and treatment for:
 - Selective destruction of RPE with ultrashort pulses (fs-ns)
 - Destruction of feeder vessels in choroidal neovascularization (CNV) related to DR
 - Other potential targets: drusen, small tumors, microaneurysms



AMD retinal pathology

Core Technology

Retinal Tracking




- **PSI's retinal tracking technology**
 - US Patents 5,767,941 and 5,943,115
 - Enables visualization, photonic analysis, and phototherapy of retinal structures with superior sensitivity and resolution
 - Absolutely required for precision laser retinal surgery




Core Technology

Adaptive Optics

- Correction of ocular aberrations by wavefront compensation & correction allows increased transverse resolution for several imaging technologies (e.g. SLO and OCT)
- For precision laser retinal surgery, AO is necessary to achieve a diffraction-limited spot at the retina



Without Compensation



With Adaptive Compensation

5 arcmin

**Cone mosaic imaged
with and without
wavefront
compensation**


**Used with permission
from Austin Roorda:
<http://www.opt.uh.edu/research/aroorda>**



Core Technology

Ultrashort Pulse Laser Damage


- Non-linear mechanisms (primarily optical breakdown) mediate retinal damage from ultrashort pulses (< 10 ns).
- Ultrashort pulse damage is therefore characterized by
 - Low energy thresholds
 - Highly localized damage with little collateral damage
- Recent engineering advances are moving ultrashort pulse lasers closer to routine use in clinical systems (OCT, refractive surgery, etc.)



Histology and OCT image of RPE disruption from 3 ps pulses.


Roach et al. JBO, in press.

System Overview



- **Dichroic beamsplitter (D)**
 - 400-860 nm $R > 95\%$
 - 870-1100 nm $T > 85\%$
- **Cold mirror (CM)**
 - 400-700 nm $R > 95\%$
 - 780-1250 nm $T > 85\%$

TAOSLO Optical Schematic



[illegible]

- **Custom PCB**
 - LUT to linearize SLO image resonant scanner
 - Synchronize and drive vertical scanner
- **3 framegrabbers (FG)**
 - TSLO
 - Flying-spot SLO
 - WS
- **RT controller**
 - Retinal tracking bandwidth > 1 kHz
 - AO closed-loop bandwidth > 10 Hz

Research Plan and Acknowledgments

- **Research Plan**

- Optical and instrumentation design complete
- System construction complete by Fall
- Human and animal subject testing in 2005

- **Acknowledgments**

- **Physical Sciences Inc.**

- R. Daniel Ferguson
- Nick Iftimia

- **Brooks AFB**

- Benjamin A. Rockwell
- Clarence Cain
- David Stolarski

- **Consulting support**

- Austin Roorda, University of Houston School of Optometry
- H. Grady Rylander III, University of Texas at Austin and Eye Institute of Austin

