

Optical Coherence Domain Reflectometry in Brain Probes


David Huang, MD, PhD

Mark Shure, PhD

Ken Baker, PhD

Ali Rezai, MD

THE CLEVELAND CLINIC





Supported by NIH/NIBIB R21 EB002718




Deep Brain Stimulation (DBS) Brain Pacemakers

FDA approved indications for
Parkinson's disease
Tremor
Dystonia

>25,000 implants to date
Reimbursed by Medicare and insurance in all states





Implantation of Deep Brain Probe




- Planning by MRI with brain atlas registration
- Stereotactic brain probe insertion guided by microelectrode recording
- Multiple insertion passes needed
- Blood vessel penetration & hemorrhagic stroke is a possible complication
- Additional intra-operative guidance mechanism needed

Preliminary Results: OCT of Rat Brain *in vitro*

- White matter tracts (axons) distinguishable from gray matter
- Greater scattering
- Higher signal
- Faster attenuation



OCDR Brain Probe Design



- 30 gauge ($D=0.3$ mm) for *in vivo* rat experiments
- 1.3/1.4 mm dual wavelength hydration measurement
- Doppler blood flow detection
- Polarization-insensitive detection

