


**Summit on Management of Radiation Dose in CT:
Toward the Sub-mSv Exam
Bethesda North Marriott, Feb. 24-25, 2011**

Innovations Required in Hardware



Willi A. Kalender, Ph.D.

**Institute of Medical Physics
University of Erlangen**

<http://www.imp.uni-erlangen.de>



Disclosures

- **WAK is a consultant to
Siemens Healthcare, Erlangen, Germany**
- **WAK is founder and shareholder of
Artemis Imaging GmbH, Erlangen, Germany**
- **WAK is a firm believer in the future of CT.**



Innovations Required in Hardware

- Definition of goals
- X-ray sources
- X-ray detectors
- Dose management
- X-ray beam collimation
- Conclusions


Estimates of effective dose in CT

Table 13: Mean effective doses for the 'Top 20 Exams' in the ten DOSE DATAMED countries

Exam type	Mean E per examination (mSv)										
	LU -	BE 00-05	DE 92-05	NO 85-95	CH 1998	FR 01-03	SE 1995	DK 1995	NL 2002	UK 90-01	Max Min
13. CT head	2.6	2.3	2.6	1.8	2.2	1.8*	2.0	1.9	1.2	2.0	2.2
14. CT neck	2.5	-	2.5	3.4	3.1	2.5*	-	1.3	-	2.4	2.6
15. CT chest	10.0	4.1	7.6	11.5	8.8	5.5*	-	11.0	5.5	7.8	2.8
16. CT spine	9.0	-	2.9	4.3	9.1	4.0*	-	5.7	3.1	4.2	3.1
17. CT abdomen	15.0	11.3	18.6	12.6	8.4	5.8*	-	14.0	10.6	9.8	3.2
18. CT pelvis	-	-	10.6	9.3	7.0	-	-	8.3	7.4	9.8	1.5
19. CT trunk	7.9	-	24.4	-	-	-	10	15.0	-	10.4	3.1
All CT	7.4	7.7	8.1	6.1	6.0	3.5*	6.0	5.9	5.3	5.4	2.3

***Natural background radiation: 3 mSv/y.
(range: 1-10 mSv/y.)***

Estimates of effective dose in CT



If dose distribution
is known
→ Organ dose
and eff. Dose E


Scan parameters
(CTDI, DLP)
are known
→ $k = E/DLP$

In general: Effective Dose $E = k \times DLP$



Multisection CT Protocols: Sex- and Age-specific Conversion Factors Used to Determine Effective Dose from Dose-Length Product¹

Paul D. Deak, PhD
Yulia Smal, MSc
Willi A. Kalender, PhD


Conversion factors CF
to estimate
effective dose E
from the dose length
product DLP
for modern scanners
using the ICRP 103
tissue weighting factors.



Dual Source CT at high pitch



- 63 y, male, 57 b.p.m.
- Pitch 3.2
- E = 0.84 mSv



European Heart Journal (2010) 31, 340–346
doi:10.1093/eurheartj/ehp470

CLINICAL RESEARCH
Imaging

Coronary computed tomography angiography with a consistent dose below 1 mSv using prospectively electrocardiogram-triggered high-pitch spiral acquisition

Stephan Achenbach^{1*}, Mohamed Marwan¹, Dieter Ropers¹, Tiziano Schepis¹, Tobias Pflederer¹, Katharina Anders², Axel Kuettner², Werner G. Daniel¹, Michael Uder², and Michael M. Leff²

¹Department of Cardiology, University of Erlangen, Ulmenweg 18, 91054 Erlangen, Germany; and ²Institute of Radiology, University of Erlangen, Germany

Received 29 May 2009; revised 1 September 2009; accepted 21 September 2009; online publish-ahead-of-print 5 November 2009



State of the Art & Goals

- Sub-mSv scanning has become a reality for a few applications already, e.g. in cardiac and pediatric CT.
- Effective dose in CT is quoted typically as 1 to 10 mSv per exam.
I here assume an average value of 5 mSv.
- To reach the goal of sub-mSv scanning in general would require a reduction by at least a factor of 5,
e.g. from 100 % = 5 mSv to 20 % = 1 mSv!





Innovations Required in Hardware

- Definition of goals
- **X-ray sources**
 - optimal spectra
 - requirements on power & filtration
- X-ray detectors
- Dose management
- X-ray collimation

Spectral Optimization for Thoracic CT

Simulations and Measurements



S: 300 x 200




M: 350 x 250

Contrast due to	Density		Iodine		Calcium	
	S	M	S	M	S	M
Size						
Optimum tube voltage	110 kV	120 kV	50 kV	60 kV	50 kV	70 kV
Change in dose at const. CNR						
120 kV → 80 kV	+ 9 %	+ 21 %	- 53 %	- 45 %	- 37 %	- 24 %

Spectral Optimization for Pediatric CT

Cadaver Measurements



Contrast due to	Iodine	Calcium
Optimum tube voltage	< 80 kV	< 80 kV
Change in dose at const. CNR		
120 kV → 80 kV	- 67 %	- 62 %

Conclusions on spectra

- Adapt the spectrum to the task and to the size of the patient to be imaged.
- There is an amazingly high potential for dose reduction by spectral optimisation.
- Higher pre-filtration helps reduce dose.
- Demands on x-ray power will increase considerably.

Dose reduction potential

	%	factor
• X-ray sources	20	0.8
• X-ray detectors		
• Dose management		
• X-ray beam collimation		


Innovations Required in Hardware

- Definition of goals
- X-ray sources
- **X-ray detectors**
- Dose management
- X-ray beam collimation
- Conclusions

Goals for detector developments in CT


- Efficiency close to 100%
for x-ray absorption and for geometry
- High spatial & temporal resolution
- Single photon counting
- Energy discrimination
- Direct conversion of x-rays
to electric signal

e.g.,
Medipix detector
CdZTe sensor
3 energy levels




Detection principles


Scintillator




Photodiode + Transistor array




Scintillator (structured)




Photodiode + Transistor array



Direct converter



Transistor array



Point Spread Functions (PSF)


Standard 64-slice Multi-Row CT detector

grid
anti-scatter
protect gaps


scintillator

reflector

optical coupling photodiode



Geometric efficiency



Flat detectors can offer 100% geometric efficiency since there are no discrete detector pixels and septa.

Absorption efficiency, %

	Attenuated spectrum (Unattenuated spectrum)					
Detector material	CsI		CdTe		GOS	
Tube voltage	80 kV	1120 kV	80 kV	1120 kV	80 kV	1120 kV
Thickness 1.00 mm	85.6	67423	99802	79900	97.8	89952
1.40 mm	92.3	79648	99903	88975	99.3	99906



Dose reduction potential

	%	factor
• X-ray sources	20	0.8
• X-ray detectors	30	0.7
• Dose management		
• X-ray beam collimation		

Innovations Required in Hardware

- Definition of goals
- X-ray sources
- X-ray detectors
- **Dose management**
 - tube current modulation (TCM)
 - automatic exposure control (AEC)
 - choice of spectra
- X-ray beam collimation
- Conclusions




Tube current modulation (TCM)

introduced in the late 1990ies
and generally available today.

The principle simply is:



Adapt the tube current and thereby the x-ray
intensity to the attenuation given for any projection
as a function of rotation angle α and position z .

In-plane or rotational TCM (α -TCM)




Attenuation for the central ray:


in a.p. direction: 43
in lateral direction: 3603



mAs Reduction by Tube Current Modulation




Conventional scan: 327 mAs




Online current modulation: 166 mAs

53% mAs reduction on average for the shoulder region
49% mAs reduction in this case

Reduction of mAs for different regions in direct comparison




mAs Reduction vs. Patient Dose Reduction




Tube Current Modulation (TCM) and Automatic Exposure Control (AEC)

Standard CT



TCM & AEC




Principle of operation

Detectors with smaller z-extent yield better performance!


Kalender WA. Computed Tomography. 2nd ed. Wiley, New York 2005



Tube Current Modulation (TCM) and Automatic Exposure Control (AEC)



Constant tube current



AEC

Resulting 3D dose distributions

**34%
exposure
reduction**


**45%
dose
reduction
(center)**



Tube current modulation (TCM)

can be applied as a function of

- rotation angle α (in-plane)
- longitudinal position z
- organ at risk, e.g. breast
- heart phase (ECG)



Conclusions on TCM and AEC

- Adapt the tube current to the size of the patient to be imaged.
- The choice of voltage and filtration, „auto-kV“, has to be a part of dose management.
- Make sure it becomes widely available and actually used.

Dose reduction potential

	%	factor
• X-ray sources	20	0.8
• X-ray detectors	30	0.7
• Dose management	30	0.7
• X-ray beam collimation		



Innovations Required in Hardware

- Definition of goals
- X-ray sources
- X-ray detectors
- Dose management
- **X-ray beam collimation**
 - dynamic z-overscanning shield
 - volume-of-interest (VOI) scanning
- Conclusions



Overlapping exposure in cone-beam CT

“step & shoot scanning”




Overscanning exposure in cone-beam CT

“spiral scanning”



- z-overscanning is associated with spiral scanning due to the need for z-interpolation.
- Wider detectors yield higher z-overscanning effects
- Up to 35% higher dose between spiral ($p=1$) and contiguous axial scans.¹

unnecessary exposure
due to z-overscanning




¹Tzedakis et al. Med Phys 2005; 20:353-59


Overcanning and Counter-Measures




Dose Profiles in Air




Dose Profiles in 32 cm CTDI Phantom



Results: Dose Saving in Air




Volume of Interest (VOI) Imaging




1. Overview (OV) scan

- complete object
- low dose
- low image quality
- VOI selection



2. Volume of Interest (VOI) scan


- increased collimation
- increased magnification
- second isocenter
- small VOI
- high local dose
- high resolution
- low noise



VOI Imaging of a breast phantom (Sim.!!)

Overview scan


0.75 mm voxels,
1.5 mGy AGD



VOI scan

left: 0.75 mm voxels,
right: 0.05 mm voxels,
1.8 mGy AGD

Cumulative dose
3.3 mGy AGD



Conclusions on beam collimation

- Limit the exposure to the volume to be imaged by dynamic collimation in z- and ϕ -direction.
- Make sure it becomes widely available and is actually used.

Dose reduction potential

	%	factor
• X-ray spectra	20	0.8
• X-ray detectors	30	0.7
• Dose management	30	0.7
• X-ray beam collimation	20	0.8
• Image reconstruction	40	0.6
Total result	80	0.2

Dose-efficient Image Reconstruction



Standard B40



SAFIRE S40

slice 1.0 mm
W = 300
C = 60

State of the art & Goals

- Sub-mSv scanning has become a reality for a few applications already.
- The goal of sub-mSv CT scanning, can be reached within this decade!
- It is not only the hardware, training & education are a necessity!
- There will be new CT applications, e.g. dedicated CT of the breast, and they will operate in the sub-mSv range!



Thank you for your attention!



willi.kalender@imp.uni-erlangen.de
www.imp.uni-erlangen.de